





AccuCore® Single Action Biopsy Needles


- Lightweight and compact – Easy to maneuver
- Sharp cutting cannula for intact core sample
- Specimen notch is fully exposed when the stylet is advanced for maximum sample
- Ultrasound enhanced tip for easier visualization under ultrasound guidance
- Depth marks on cannula for precise positioning
- 14 ga has Trocar point for ease of insertion in breast tissue



AccuCore® Single Action Biopsy Needle



Standard Point



Trocar Point
14 ga only

AccuCore® Single Action Biopsy Needle – Catalog Codes

Code	Gauge	Length (CM)	Length (In)	Tip	Coaxial w/HiLiter® Ultrasound Tip
640914	14	9	3½"	Trocar	601409
641514	14	15	5⅞"	Trocar	601415
620916	16	9	3½"	Standard	601609
621516	16	15	5⅞"	Standard	601615
620918	18	9	3½"	Standard	601809
621518	18	15	5⅞"	Standard	601815
620920	20	9	3½"	Standard	602009
621520	20	15	5⅞"	Standard	602015

Patent No. 4,582,061, 5,490,521, and 5,820,554.

INRAD Inc., 4375 Donker Ct. SE, Grand Rapids, MI 49512, Tel: 800-558-4647, 616-301-7800, Fax: 616-301-7799, E-mail: inrad-inc@inrad-inc.com

06369C

For More Information Call Toll Free:
1-800-558-4647



4375 Donker Ct. SE
Grand Rapids, MI 49512

Toll Free: 800-558-4647
Tel: 616-301-7800
Fax: 616-301-7799
Email: inrad-inc@inrad-inc.com
www.inrad-inc.com

Complete information on INRAD devices, including indications, contraindications and cautions, is contained in the product package labeling and should be consulted before use.

Restricted Device: Federal law (USA) restricts the sale, distribution or use of this device to, by or on the lawful order of a physician. This product is designed as a disposable product for single patient use only.

® INRAD is a registered trademark of INRAD, Inc.
™ Ghiatas is a trademark of INRAD, Inc.
® SampleMASTER is a registered trademark of MEDTRONIC, Inc.
® HiLiter is a registered trademark of MEDTRONIC, Inc.
® SureStop is a registered trademark of MEDTRONIC, Inc.
® AccuMark is a registered trademark of MEDTRONIC, Inc.
™ Guard Guide is a trademark of INRAD, Inc.
™ Procytin is a trademark of INRAD, Inc.
® AccuCore is a registered trademark of INRAD, Inc.
® PreSet is a registered trademark of INRAD, Inc.
® Express is a registered trademark of INRAD, Inc.
® UltraClip is a registered trademark of INRAD, Inc.
® ViewPoint is a registered trademark of INRAD, Inc.
™ SentiLoc is a trademark of INRAD, Inc.
™ UltraWire is a trademark of INRAD, Inc.
™ Make Your Mark is a trademark of INRAD, Inc.
® AccuPlace is a registered trademark of MEDTRONIC, Inc.
© Copyright 2002-2006 INRAD, Inc.