




# Chiba Type Biopsy Needles

## All INRAD® Chiba Type Biopsy Needles:


- Crystal clear hub
- Color coded according to gauge
- Beveled outer cannula provides excellent collections
- Available with AccuMark® Sheath for easy depth marking

## SampleMASTER® Series has these additional features:

- SureStop® depth marker that stays where you put it
- Locking hub for ease of handling, no stylet and cannula separation until you want it
- Depth marks etched on cannula
- Palpable alignment markings on hub for easy orientation
- Optional HiLiter® Ultrasound Enhancement, extremely visible



**Chiba Needle**



**AccuMark SHEATH**

**Chiba Needle Point Detail**



**SampleMASTER Series**

## Chiba Type Biopsy Needle – Catalog Codes

Description	Standard	With AccuMark Sheath	SampleMASTER	SampleMASTER with HiLiter Ultrasound Enhancement
18 ga x 2 3/8" (6.0 cm)	NA	NA	NA	53218
18 ga x 5 7/8" (14.9 cm)	54018	45018	54818	53818
18 ga x 7 7/8" (20.0 cm)	54118	45118	NA	NA
20 ga x 2 3/8" (6.0 cm)	NA	NA	NA	53220
20 ga x 3 7/8" (9.8 cm)	NA	NA	54720	53720
20 ga x 5 7/8" (14.9 cm)	54020	45020	54820	53820
20 ga x 7 7/8" (20.0 cm)	54120	45120	NA	NA
22 ga x 2 3/8" (6.0 cm)	NA	NA	NA	53222
22 ga x 3 7/8" (9.8 cm)	54222	45222	54722	53722
22 ga x 5 7/8" (14.9 cm)	54022	45022	54822	53822
22 ga x 7 7/8" (20.0 cm)	54122	45122	54922	53922
Packaged	20 per carton	20 per carton	20 per carton	10 per carton

US Patent No. 4582061, 5217438, 5336191, 5490521, 5820554, 5466225, 5571091, and others pending.

For More Information Call Toll Free:  
**1-800-558-4647**



4375 Donker Ct. SE  
Grand Rapids, MI 49512

Toll Free: 800-558-4647  
Tel: 616-301-7800  
Fax: 616-301-7799  
Email: [inrad-inc@inrad-inc.com](mailto:inrad-inc@inrad-inc.com)  
[www.inrad-inc.com](http://www.inrad-inc.com)

Complete information on INRAD devices where applicable, including indications, contraindications and cautions, is contained in the product package labeling and should be consulted before use.

Restricted Device: Federal law (USA) restricts the sale, distribution or use of this device to, by or on the lawful order of a physician. This product is designed as a disposable product for single patient use only.

® INRAD is a registered trademark of INRAD, Inc.  
™ Ghiatas is a trademark of INRAD, Inc.  
® SampleMASTER is a registered trademark of MEDTRONIC, Inc.  
® HiLiter is a registered trademark of MEDTRONIC, Inc.  
® SureStop is a registered trademark of MEDTRONIC, Inc.  
® AccuMark is a registered trademark of MEDTRONIC, Inc.  
™ Guard Guide is a trademark of INRAD, Inc.  
™ Procytin is a trademark of INRAD, Inc.  
® AccuCore is a registered trademark of INRAD, Inc.  
® PreSet is a registered trademark of INRAD, Inc.  
® Express is a registered trademark of INRAD, Inc.  
® UltraClip is a registered trademark of INRAD, Inc.  
® ViewPoint is a registered trademark of INRAD, Inc.  
™ SentiLoc is a trademark of INRAD, Inc.  
™ UltraWire is a trademark of INRAD, Inc.  
™ Make Your Mark is a trademark of INRAD, Inc.  
® AccuPlace is a registered trademark of MEDTRONIC, Inc.  
© Copyright 2002-2006 INRAD, Inc.