




## Coaxial Biopsy Needles

- Allows multiple sample retrieval from single entry site
- Choice of Franseen or Chiba biopsy needle
- Ultrasound enhanced introducer and biopsy needle
- SureStop® depth marker that stays where you put it
- Depth marks on introducer
- Depth marks on biopsy needles for precise placement
- 19 ga x 10 cm introducer
- 22 ga x 14.9 cm biopsy needle



**Introducer Needle**



### Coaxial Biopsy Needle – Catalog Codes

Description	Catalog Codes
<b>Franseen Biopsy Needle</b> 22 ga x 5 <sup>7</sup> / <sub>8</sub> " (14.9 cm) 19 ga x 4" (10.0 cm) Introducer Needle	54622
<b>Chiba Biopsy Needle</b> 22 ga x 5 <sup>7</sup> / <sub>8</sub> " (14.9 cm) 19 ga x 4" (10.0 cm) Introducer Needle	54522

Packaged 10 per carton

Patent No. 5217438

INRAD Inc., 4375 Donker Ct. SE, Grand Rapids, MI 49512, Tel: 800-558-4647, 616-301-7800, Fax: 616-301-7799, E-mail: inrad-inc@inrad-inc.com

06220B

For More Information Call Toll Free:  
**1-800-558-4647**



4375 Donker Ct. SE  
Grand Rapids, MI 49512

Toll Free: 800-558-4647  
Tel: 616-301-7800  
Fax: 616-301-7799  
Email: [inrad-inc@inrad-inc.com](mailto:inrad-inc@inrad-inc.com)  
[www.inrad-inc.com](http://www.inrad-inc.com)

Complete information on INRAD devices, including indications, contraindications and cautions, is contained in the product package labeling and should be consulted before use.

Restricted Device: Federal law (USA) restricts the sale, distribution or use of this device to, by or on the lawful order of a physician. This product is designed as a disposable product for single patient use only.

® INRAD is a registered trademark of INRAD, Inc.  
™ Ghiatas is a trademark of INRAD, Inc.  
® SampleMASTER is a registered trademark of MEDTRONIC, Inc.  
® HiLiter is a registered trademark of MEDTRONIC, Inc.  
® SureStop is a registered trademark of MEDTRONIC, Inc.  
® AccuMark is a registered trademark of MEDTRONIC, Inc.  
™ Guard Guide is a trademark of INRAD, Inc.  
™ Procytin is a trademark of INRAD, Inc.  
® AccuCore is a registered trademark of INRAD, Inc.  
® PreSet is a registered trademark of INRAD, Inc.  
® Express is a registered trademark of INRAD, Inc.  
® UltraClip is a registered trademark of INRAD, Inc.  
® ViewPoint is a registered trademark of INRAD, Inc.  
™ SentiLoc is a trademark of INRAD, Inc.  
™ UltraWire is a trademark of INRAD, Inc.  
™ Make Your Mark is a trademark of INRAD, Inc.  
® AccuPlace is a registered trademark of MEDTRONIC, Inc.  
© Copyright 2002-2006 INRAD, Inc.