



Franseen Type Biopsy Needles

All INRAD® Franseen Type Biopsy Needles:


- Crystal clear hub
- Color coded according to gauge
- Notched crown shape aids in cutting and retrieving of cells and small cores of tissue
- Available with AccuMark® Sheath for easy depth marking

SampleMASTER® Series has these additional features:


- SureStop® depth marker that stays where you put it
- Locking hub for ease of handling, no stylet and cannula separation until you want it
- Depth marks etched on cannula
- Palpable alignment markings on hub for easy orientation
- HiLiter® Ultrasound Enhancement, extremely visible



Franseen Needle



Franseen Needle Point Detail



AccuMark SHEATH



SampleMASTER Series

Franseen Type Biopsy Needle – Catalog Codes

Description	Standard	With AccuMark Sheath	SampleMASTER with HiLiter Ultrasound Enhancement
18 ga x 7 7/8" (20.0 cm)	42118	NA	NA
20 ga x 5 7/8" (14.9 cm)	42020	44020	41820
20 ga x 7 7/8" (20.0 cm)	42120	44120	NA
22 ga x 5 7/8" (14.9 cm)	NA	44022	41822
22 ga x 7 7/8" (20.0 cm)	42122	NA	NA
Packaged	20 per carton	20 per carton	10 per carton

US Patent No. 4582061, 5217438, 5336191, 5490521, 5820554, 5466225, 5571091, and others pending.

For More Information Call Toll Free:
1-800-558-4647



4375 Donker Ct. SE
Grand Rapids, MI 49512

Toll Free: 800-558-4647
Tel: 616-301-7800
Fax: 616-301-7799
Email: inrad-inc@inrad-inc.com
www.inrad-inc.com

Complete information on INRAD devices where applicable, including indications, contraindications and cautions, is contained in the product package labeling and should be consulted before use.

Restricted Device: Federal law (USA) restricts the sale, distribution or use of this device to, by or on the lawful order of a physician. This product is designed as a disposable product for single patient use only.

® INRAD is a registered trademark of INRAD, Inc.
™ Ghiatas is a trademark of INRAD, Inc.
® SampleMASTER is a registered trademark of MEDTRONIC, Inc.
® HiLiter is a registered trademark of MEDTRONIC, Inc.
® SureStop is a registered trademark of MEDTRONIC, Inc.
® AccuMark is a registered trademark of MEDTRONIC, Inc.
™ Guard Guide is a trademark of INRAD, Inc.
™ Procytin is a trademark of INRAD, Inc.
® AccuCore is a registered trademark of INRAD, Inc.
® PreSet is a registered trademark of INRAD, Inc.
® Express is a registered trademark of INRAD, Inc.
® UltraClip is a registered trademark of INRAD, Inc.
® ViewPoint is a registered trademark of INRAD, Inc.
™ SentiLoc is a trademark of INRAD, Inc.
™ UltraWire is a trademark of INRAD, Inc.
™ Make Your Mark is a trademark of INRAD, Inc.
® AccuPlace is a registered trademark of MEDTRONIC, Inc.
© Copyright 2002-2006 INRAD, Inc.