

# Select Core™

## Variable Throw Biopsy Device

*Designed for ease of  
use and flexibility in  
soft tissue biopsy*



**Fully  
Adjustable  
Throw**



**Manual or  
Fully Automatic**



*...helping to reduce the uncertainty in diagnosis.*



*Variable Throw Indicator*



*Thumb Wheel for Throw Adjustment*

## Variable Throw

The SelectCore™ features the ability to set the precise throw needed for the specimen size desired or to avoid critical structures during the biopsy. After arming, simply rotate the thumb wheel at the front of the device to any distance between 10 and 25mm.



## Two Firing Modes

For flexibility, the SelectCore™ offers a choice of either manual or fully automatic modes. Manual mode deploys the stylet only and is ideal for denser tissue or when there is a desire to verify notch position before tissue capture. Automatic mode fires both the stylet and cannula in sequence for tissue capture.

## Additional Features

- Ergonomic design allows for ease of use along with comfort and control
- Easy to grip exterior and ring-shaped pull lever handle for straightforward device arming
- Ultrasound enhanced tip for increased visualization under ultrasound guidance
- Depth marks etched on cannula for precise positioning
- Color coded for quick gauge size identification
- Matching coaxial introducer needle featuring HiLiter® ultrasound enhanced tip available



Ordering Information				
SelectCore™ Variable Throw Biopsy Device				
Code	Gauge	Length (CM)	Length (In)	Coaxial w/ HiLiter® Tip
991014	14	10	4"	601409
991514	14	15	5 7/8"	691514
991018	18	10	4"	601809
991518	18	15	5 7/8"	691518
992018	18	20	7 7/8"	692018
991020	20	10	4"	602009
991520	20	15	5 7/8"	691520
992020	20	20	7 7/8"	692020

**For More Information Call Toll Free:  
1-800-558-4647**



4375 Donker Ct. SE  
Grand Rapids, MI 49512

Toll Free: 800-558-4647  
Tel: 616-301-7800  
Fax: 616-301-7799  
Email: [inrad-inc@inrad-inc.com](mailto:inrad-inc@inrad-inc.com)  
[www.inrad-inc.com](http://www.inrad-inc.com)

Complete information on INRAD devices, including indications, contraindications and cautions, is contained in the product package labeling and should be consulted before use.

Restricted Device: Federal law (USA) restricts the sale, distribution or use of this device to, by or on the lawful order of a physician. This product is designed as a disposable product for single patient use only.

©INRAD is a registered trademark of INRAD, Inc.  
©HiLiter is a registered trademark of Medtronic, Inc.  
™SelectCore is a trademark of INRAD, Inc.  
© Copyright 2009 INRAD, Inc.  
Patents Pending  
01902A